

Put your trust in the proven leader of regulator design and development with 50 plus years of customer satisfaction! The unique design of the TRG Heliox Regulator will dramatically improve treatment of asthma emergencies! The Diameter Index Safety System (DISS) Fitting allows for the operation of a nebulizer simultaneously with TRG's Heliox Regulator for improved drug delivery and more successful asthma emergency rescues. Proudly made in the U.S.A.

## Heliox



### SPECIFICATIONS

Precise Flow Settings from 0-20 lpm  
 Helium/Oxygen Ratio of 80:20 or 70:30  
 Available With or Without a Diameter Index Safety System (DISS 1180) Fitting  
 CGA 890 or CGA 280  
 Ask about complete kit purchases



Part Num.	Discriptions
2210-20	70:30 Helium/Oxygen ratio, CGA 890, yoke style, 0-20 LPM
2211-20	80:20 Helium/Oxygen ratio, CGA 890, yoke style, 0-20 LPM
2212-20	70:30 Helium/Oxygen ratio, CGA 280, nut & nipple style, 0-20 LPM
2213-20	80:20 Helium/Oxygen ratio, CGA 280, nut & nipple style, 0-20 LPM
2214-20	70:30 Helium/Oxygen ratio, CGA 890, yoke style with DISS fitting, 0-20 LPM
2215-20	80:20 Helium/Oxygen ratio, CGA 890, yoke style with DISS fitting, 0-20 LPM
2216-20	70:30 Helium/Oxygen ratio, CGA 280, nut & nipple style with DISS fitting, 0-20 LPM
2217-20	80:20 Helium/Oxygen ratio, CGA 280 nut & nipple style with DISS fitting, 0-20 LPM

TRG is a low-cost, high-volume manufacturer of medical regulators, oxygen conserving devices, portable liquid oxygen systems, the exclusive distributor of Piper Impact aluminum medical cylinders and more. TRG's products are used all over the world by some of the largest and well respected homecare and hospital providers. As a manufacturer, TRG is able to meet specific customer needs while maintaining competitive price structures, stringent quality systems, and a knowledgeable customer care staff. TRG goes above and beyond for you...

*every breath of the way...*