

The Escort Nebulizer is a stationary compressor style device for individuals with Asthma, COPD, or other respiratory disorders. Trusted to deliver concise and reliable, piston-driven treatment of all popular solution medications.

- Nebulizer Capacity: 6 ml (cc)
- Output aerosol: 1.91 mL
- Particulate size: 0.5 - 5microns
- Output rate: approximately 0.16 ml/min
- Compressor output flow rate: 10 ± 2 LPM
- Operating voltage: 120 VAC @ 60 Hz
- Power consumption: 1.4A/68W
- Compressor noise level: less than 58 dB
- Approximate operating pressure: 18 psi @ 6 LPM
- Nebulizer efficiency: 63.7%
- Rate Efficiency Index: 0.102
- Double-insulated Class II
- 510K approved
- Limited 2 year warranty

Escort Nebulizer



6 ml Capacity and Very High Efficiency!



Convenient Power Cord Storage Cavity



Accessories Included/Replacements Avail.



Comfortable Face Forming Mask



Easy to Clean Mouth Piece and Adaptors

TRG is a low-cost, high-volume manufacturer of medical regulators, oxygen conserving devices, portable liquid oxygen systems, the exclusive distributor of Piper Impact aluminum medical cylinders and more. TRG's products are used all over the world by some of the largest and well respected homecare and hospital providers. As a manufacturer, TRG is able to meet specific customer needs while maintaining competitive price structures, stringent quality systems, and a knowledgeable customer care staff. TRG goes above and beyond for you...

every breath of the way...